





Dale NasoGastric Tube Holder

APPLICATION INSTRUCTIONS

APPLICATION INSTRUCTIONS

1. After the tube is in place, make certain that the nose and the NG tube are both clean and dry.
2. Cleanse the nose and area around it with a clean wash cloth that has been moistened with warm clean tap water(*NO SOAP*). Towel dry the area immediately after cleaning with a clean dry towel. (*DO NOT APPLY ANY ALCOHOL TO THIS APPLICATION AREA*) Alcohol will dry the tissue and may later cause micro-cracking of the skin.(Figure A)
3. We also recommend cleaning the exterior of the NG tube below the nare area to remove and residual lubrication gel prior to wrapping with the NG Tube Holder.
4. Peel back the NG Tube Holder from the liner starting at the top of the nose pad.
5. Position the nosepad APEX(^) at a slight angle towards the nare holding the NG tube on the patient's nose, at approximately ¼" to ½" below the tip of the nose. This will help allow the tube to hang freely.(Figure B)
6. Choose the tab closest to the tube and begin to wrap it in a spiraling fashion. Try and leave a small window at the beginning of the wrap to expose the cm marking of the NG tube.(Figure C)
7. Wrap the second remaining tab using the same spiral technique.(Figure D)
8. You may mark the date the NGTH was applied with a non-toxic pen or marker.
9. Change the NG Tube Holder after 3 days, or sooner if necessary.



REMOVAL INSTRUCTIONS

1. To remove the NG tube holder from the patient, simply pull back on the non-adherent tips on the end of each tab to unwrap them from the tube.
2. Continue to gently pull the tube holder in an upward fashion to remove from the patient's nose.

* This is a single patient use disposable product

ADDITIONAL BENEFITS

- Ideal for holding most nasogastric tubes, including nasoenteric, Salem Sump™, and levin stomach tubes.
- Skin-friendly adhesive stays in place for up to 3 days, yet is easy to remove.
- Non-adherent plastic tips on both tabs make tube removal safer and easier when donning gloves.

Product #	Description	Quantity
160	Large / Adult	50/box
161	Small / Young Adult / Pediatric / Geriatric	50/box

CLINICAL STUDIES

The Dale NasoGastric Tube Holder was used in a recent study published in MEDSURG Nursing Magazine. To view the article please visit dalemed.com/ng-study.



SCAN FOR
MORE INFO

