

Hold-n-Place® Multipurpose Tube Holder



The Dale Hold-n-Place Multipurpose Tube Holder has a skin friendly, breathable, adhesive base that secures a wide variety of catheters and tubes from 12fr up to 40fr, including: Jackson-Pratt, PEG-tubes, suprapubic tubes, IV lines, chest tubes, Foley catheters, drainage tubing, and multi-port catheters.

Simplified Closure System

The Holder features a unique hook-and-loop closure that secures tubing without restricting movement or irritating the skin. Its simple design speeds application and adjustment.

The center section includes optional adhesive:

- + **With liner left on** – provides secure, adhesive-free retention.
- + **With liner removed** – offers adhesive locking around the tube.

Compatible with all tubing types, including Teflon™, silver, and silicone-coated catheters.

Key Features

- ◆ Compatible with a wide range of catheters and tubes (12Fr-40Fr), including Teflon™, silver, and silver-coated options for standardized use.
- ◆ Skin-friendly, breathable adhesive base allows secure wear for up to 7 days.
- ◆ Easy-to-apply base with unique closure for faster application and simple adjustments.
- ◆ Offers choice of adhesive or non-adhesive tube retention.
- ◆ DEHP-free and not made with natural rubber latex.

Product #	Description	Quantity
150	Hold-n-Place Multipurpose Tube Holder Fits catheters and tubes 12fr-40fr	50/box

HCPCS Codes: Adhesive Patch - A4333



Hold-n-Place® Multipurpose Tube Holder

Features	Clinical Advantages
Base	■ High breathability and non-absorbent surface for extended use (up to 7 days). ■ Rounded, contoured edges, and high skin adhesion strength prevent premature lift-off. ■ Sized to fit most adults and pediatrics. ■ Writeable surface.
Low Profile	■ Designed to reduce the likelihood of the tube being caught on bed or clothing. ■ No floppy closure, thereby reducing tissue irritation. ■ Helps prevent tube dislodgement.
Exclusive Closure System	■ Secures a variety of catheters and tubes, sizes 12fr-40fr. ■ No fear of choosing the incorrect size, one product fits most, reducing inventory. ■ Mating hook and loop locks the closure in place. ■ Catheter/tube can be adhesively or non-adhesively secured into position. ■ Reduces the risk of inadvertent dislodgement. ■ Stabilizes catheter/tube to prevent pistoning, ■ Re-sealable adhesive tab allows for repositioning of catheter/tubes several times. ■ Simplifies catheter/tube care.
No Hard Plastic Parts	■ Reduces likelihood of skin irritation and skin breakdown. ■ No hard plastic components to dig into skin.
Not made with natural rubber latex	■ Eliminates patient and caregiver latex reactions.

Learn More

