

# Tracheostomy Tube Holders



Dale Tracheostomy Tube Holders eliminate the frustrations associated with twill ties and other holders while minimizing secondary complications at the same time.

## Key Features

- ◆ Moisture wicking neckband lining reduces the risk of skin breakdown by minimizing the amount of moisture absorbed
- ◆ Adjustable two piece neckband
- ◆ PediPrints® series offers kids fun, comforting prints†
- ◆ Fastener tabs quickly and easily attach to any size tracheostomy tube
- ◆ Stretch material allows for cough reflex, accommodates edema, and ensures a snug fit without compromising security‡
- ◆ Not made with natural rubber latex

† On the 241 and 242 tube holder

‡ On the 240 Blue tube holder

Product #	Description	Quantity
<b>240</b>	<b>240 Blue®</b> 2.54cm(1") wide neckband, fits up to 49.5cm(19.5") neck	10/box
<b>241</b>	<b>241 PediPrints®</b> 2.54cm(1") wide neckband, fits up to 46cm(18") neck	10/box
<b>242</b>	<b>242 PediPrints®</b> 1.91cm(0.75") wide neckband, fits up to 23cm(9") neck	10/box
<b>244</b>	<b>244 Blue®</b> 3.81cm(1.5") wide neckband, fits up to 71cm(28") neck	10/box

HCPCS Code A7526 • Patents [dalemed.com/patents](http://dalemed.com/patents)



# Tracheostomy Tube Holders

Features	Clinical Advantages
<b>Quick-Fastening Tabs</b>	<ul style="list-style-type: none"><li>▪ Easy to thread through flange of any size tracheostomy plate</li><li>▪ Helps prevent movement of the trach plate</li><li>▪ Patented tapered fastener tab (244 Blue® only)</li><li>▪ Securely holds trach tube in place</li><li>▪ Allows for easy change of tube holder and trach dressings</li><li>▪ Will hold even if tabs are wet</li></ul>
<b>Moisture Wicking Lining</b>	<ul style="list-style-type: none"><li>▪ Gentle, even on delicate skin</li><li>▪ Helps keep neck dry, prevents skin excoriation on neck</li></ul>
<b>Adjustable Neckband</b>	<ul style="list-style-type: none"><li>▪ Adjusts to fit most patients, from pediatric to adult</li><li>▪ Width reduces rolling or folding</li><li>▪ Allows quick and easy adjustment of neckband</li><li>▪ Comfortable, padded neckband won't dig into skin</li></ul>
<b>Stretch Material on Neckband (240 Blue® only)</b>	<ul style="list-style-type: none"><li>▪ Accommodates edema and cough reflex</li><li>▪ Accommodates weight of ventilator circuitry</li><li>▪ Maintains a snug, comfortable fit</li></ul>
<b>PediPrints®</b>	<ul style="list-style-type: none"><li>▪ Comforting patterns for patient and caregiver</li></ul>
<b>Not Made with Natural Rubber Latex</b>	<ul style="list-style-type: none"><li>▪ Eliminates patient and caregiver latex reactions</li></ul>

Learn More

