

Nasal Dressing Holder



The Dale Nasal Dressing Holder replaces adhesive tape to hold nasal dressings in place following rhinoplasty, septoplasty, sinus surgery, facial trauma, and nose bleeds.

Key Features

- ◆ Non-adhesive holder eliminates irritation and excoriation of facial skin.
- ◆ Plastic center section holds dressings in place without absorbing secretions.
- ◆ Hook tabs can be adjusted to accommodate most patients.
- ◆ Soft ear loops simplify dressing changes.
- ◆ Saves nursing time.
- ◆ Improves patient comfort.
- ◆ One size fits most patients.
- ◆ Not made with natural rubber latex.

Product #	Description	Quantity
600	Nasal Dressing Holder One size fits most	10/box



Nasal Dressing Holder

Features	Clinical Advantages
Non-absorbent Plastic Center Section	<ul style="list-style-type: none">■ Securely holds drainage “drip” pads and packing following nasal surgeries.■ Wipeable plastic center section holds dressings in place.■ Eliminates skin breakdown and allergies related to tape.■ Enhances recovery by increasing patient confidence and sense of security with mobility.■ Soft moisture repellent material adds to comfort.
Soft, Stretchy, Earloops	<ul style="list-style-type: none">■ Eliminates irritation to the anchoring point behind the ears.■ Stretches to accommodate edema.■ Allows for easy dressing changes.
Hook Tabs	<ul style="list-style-type: none">■ Easy to reposition hook to individualize size for each patient.
Easy to Use	<ul style="list-style-type: none">■ Reduces nursing time by reducing time to change dressing.■ Allows patient to change holder without fear of pain associated with adhesive tape.
Not Made with Natural Rubber Latex	<ul style="list-style-type: none">■ Eliminates patient and caregiver latex reactions.

Learn More

