


Nasal Dressing Holder



The Dale Nasal Dressing Holder replaces adhesive tape to hold nasal dressings in place following rhinoplasty, septoplasty, sinus surgery, facial trauma, and nose bleeds.

Key Features

- ◆ Non-adhesive holder eliminates irritation and excoriation of facial skin.
- ◆ Plastic center section holds dressings in place without absorbing secretions.
- ◆ Hook tabs can be adjusted to accommodate most patients.
- ◆ Soft ear loops simplify dressing changes.
- ◆ Saves nursing time.
- ◆ Improves patient comfort.
- ◆ One size fits most patients.
- ◆ Not made with natural rubber latex.

Product #	Description	Quantity
600	Nasal Dressing Holder One size fits most	10/box



Nasal Dressing Holder

Features	Clinical Advantages
Non-absorbent Plastic Center Section	<ul style="list-style-type: none">Securely holds drainage “drip” pads and packing following nasal surgeries.Wipeable plastic center section holds dressings in place.Eliminates skin breakdown and allergies related to tape.Enhances recovery by increasing patient confidence and sense of security with mobility.Soft moisture repellent material adds to comfort.
Soft, Stretchy, Earloops	<ul style="list-style-type: none">Eliminates irritation to the anchoring point behind the ears.Stretches to accommodate edema.Allows for easy dressing changes.
Hook Tabs	<ul style="list-style-type: none">Easy to reposition hook to individualize size for each patient.
Easy to Use	<ul style="list-style-type: none">Reduces nursing time by reducing time to change dressing.Allows patient to change holder without fear of pain associated with adhesive tape.
Not Made with Natural Rubber Latex	<ul style="list-style-type: none">Eliminates patient and caregiver latex reactions.

Learn More

